

DATA SHEET

Rev.ver: 01

1. Identification

Product Name: Dulbecco's Phosphate Buffered Saline 10X w/o Phenol red, Calcium and Magnesium

Catalog.no: CCBL004

2. Description

Balanced Salt Solution (Physiological Solution) have the properties of buffer solution (regulate pH), normal saline (maintain osmotic pressure) and culture medium (provide nutrition). It can meet the basic needs of survival and metabolism of tissues, organs or cells in vitro. A small amount of phenolic red was added to some equilibrium salt solutions to indicate the pH change of the solution. Dulbecco's Phosphate Buffered Saline (DPBS) is one of the most widely used equilibrium salt solutions in biochemistry. The main components are NaCl, KCl, KH_2PO_4 and Na_2HPO_4 . DPBS was mainly used for embryological research.

3. Storage

Store below 30°C . Use before expiry date on the label.

4. Composition

Ingredients	mg/L
Potassium Chloride	2000.00
Disodium hydrogen phosphate	11500.00
Potassium dihydrogen orthophosphate	2000.00
Sodium chloride	80000.00

5. Technical Support

For technical assistance or more information, please contact: BioArtha Labs Pvt Ltd or **visit:** www.bioarthalabs.com, **Email:** info@bioarthalabs.com.

Disclaimer: This data sheet is provided for informational purposes only. The performance and stability of Dulbecco's Phosphate Buffered Saline 10X w/o Phenol red, Calcium and Magnesium may vary depending on experimental conditions. It is recommended to optimize reaction conditions for specific applications